
Silica gel 60-120 mesh (125 - 250 μm)

GRM7477

Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white fine powder or crystals
Solubility	:	Insoluble in water
pH (10% aqueous suspension at 25°C)	:	6.50 - 7.50
Particle size 1	:	min. 80.00% passing through 60 mesh (250 μm)
Particle size 2	:	max. 20.00% passing through 120 mesh (125 μm)
Particle size	:	60 - 120 mesh (125 - 250 μm)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg