
Silica gel 60 - 200 mesh (74 - 250 μm)

GRM7478

Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white fine powder or crystals
Solubility	:	Insoluble in water
pH (10% aqueous suspension at 25°C)	:	6.50 - 7.50
Loss on drying (at 150°C, 4 hr)	:	\leq 2.00%
Particle size 1	:	min. 80.00% passing through 60 mesh (250 μm)
Particle size 2	:	max. 20.00% passing through 200 mesh (74 μm)
Particle size	:	60 - 200 mesh (74 - 250 μm)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg