

---

## Silica gel 100-200 mesh (74 - 149 $\mu\text{m}$ )

**GRM7479**

---

### Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	Colourless or white fine powder or crystals
Solubility	:	Insoluble in water
Loss on drying (at 150°C, 4 hr)	:	$\leq$ 2.00%
Particle size 1	:	min. 80.00% passing through 100 mesh (149 $\mu\text{m}$ )
Particle size 2	:	max. 20.00% passing through 200 mesh (74 $\mu\text{m}$ )
Particle size	:	100 - 200 mesh (74 -149 $\mu\text{m}$ )

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg