
Silica gel 200 - 400 Mesh (37 - 74 μ m)

GRM7481

Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white fine powder or crystals
Solubility	:	Insoluble in water
pH (10% aqueous suspension at 25°C)	:	6.00 - 7.00
Loss on drying (at 105°C, 2 hr)	:	\leq 4.00%
Particle size 1	:	min. 80.00% passing through 200 mesh (74 μ m)
Particle size 2	:	max. 20.00% passing through 400 mesh (37 μ m)
Particle size	:	200 - 400 mesh (37 -74 μ m)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg