
Silica gel 230 - 400 mesh (37 - 63 μm),

GRM7484

For flash chromatography

Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white fine powder or crystals
Solubility	:	Insoluble in water
pH (10% aqueous suspension at 25°C)	:	6.50 - 7.50
Loss on drying (at 105°C, 2 hr)	:	\leq 7.00%
Particle size 1	:	min. 80.0% passing through 230 mesh (63 μm)
Particle size 2	:	max. 20.00% passing through 400 mesh (37 μm)
Particle size	:	230 - 400 mesh (37 - 63 μm)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg