

Sodium butyrate, Hi-LR[®]

GRM7497


Product Identifier

CAS No.	:	156-54-7
EC No.	:	205-857-6
Molecular Formula	:	C ₄ H ₇ NaO ₂
Molecular Weight	:	110.09
HS Code	:	2915 60 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Water (K.F.)	:	<= 5.00%
Assay (GC/NT, on anhydrous basis)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	ET6400000
WGK	:	1