

---

## di-Sodium hydrogen citrate sesquihydrate

**GRM7503**

### Product Identifier

---

CAS No.	:	6132-05-4
EC No.	:	205-623-3
Molecular Formula	:	$C_6H_6Na_2O_7 \cdot 1.5H_2O$
Molecular Weight	:	263.11
HS Code	:	2918 15 40
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colourless to off-white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	4.60 - 5.20
Assay (NT/NaOH Titration)	:	98.00 - 104.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
WGK	:	3