

---

## Succinic acid, Hi-LR™

**GRM7509**

### Product Identifier

---

CAS No.	:	110-15-6
EC No.	:	203-740-4
Molecular Formula	:	C <sub>4</sub> H <sub>6</sub> O <sub>4</sub>
Molecular Weight	:	118.09
HS Code	:	2917 19 30
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colourless to white crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	184 - 191°C
Assay (NaOH Titration)	:	min. 98.00%

### Safety Information

---

Hazard Pictogram(s) :



Signal Word	:	Danger
Hazard Statement(s)	:	H318
Precautionary Statement(s)	:	P280- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WM4900000
WGK	:	1