

GRM7514 Sodium pentachlorophenate**Product Number Packing**

GRM7514 : 500G

Product Information

Product Code : GRM7514
Product Name : Sodium pentachlorophenate
Synonym : Pentachlorophenol sodium salt
Molecular Formula : C_6Cl_5NaO
Molecular Weight : 288.32
CAS No. : 131-52-2
EC No. : 205-025-2
Shelf Life : 4 years

Technical Specification

Appearance : Buff crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
FTIR (KBr disc) : Matches with the standard pattern
Water (K.F.) : $\leq 15.00\%$
Assay (NT) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H301-H311-H315-H319-H330-H335-H351-H410
Precautionary Statement(s) : P260-P273-P280-P284-P305+P351+P338-P310
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 24/25-26-36/37/38-40-50/53
S-Phrase(s) : 22-28-36/37-45-52-60-61
WGK : 3
RTECS : SM6490000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 2567
Class : 6.1
Packaging Group : 2
Marine Pollutant : Yes



Product Information

ADR/RID : 2567 6.1/PG 2
IMDG : 2567 6.1/PG 2
IATA : 2567 6.1/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein