

tetra-sodium pyrophosphate anhydrous

GRM7515


Product Identifier

CAS No.	:	7722-88-5
EC No.	:	231-767-1
Molecular Formula	:	Na ₄ P ₂ O ₇
Molecular Weight	:	265.90
HS Code	:	2835 39 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	9.50 - 10.80
Assay (HCl Titration)	:	97.50 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H318
Precautionary Statement(s)	:	P280- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UX7350000
WGK	:	1