

GRM7517**Sodium salicylate****Product Number****Packing**

GRM7517 : 500G

Product Information

Product Code : GRM7517
Product Name : Sodium salicylate
Synonym : 2-Hydroxybenzoic acid sodium salt; salicylic acid sodium salt
Molecular Formula : $C_7H_5NaO_3$
Molecular Weight : 160.10
CAS No. : 54-21-7
EC No. : 200-198-0
HS Code : 2918 21 00
Shelf Life : 4 years

Technical Specification

Appearance : Colourless or white to off-white crystals or powder or flakes
Solubility : 33.3 mg soluble in 1 mL of water
Assay (NT) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302-H319
Precautionary Statement(s) : P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22-36
S-Phrase(s) : 26
WGK : 1
RTECS : VO5075000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods