

RM7519 Sodium starch glycolate, Hi-LR™

Product Number Packing

RM7519 : 500G
RM7519 : 5KG

Product Information

Product Code : RM7519
Product Name : Sodium starch glycolate, Hi-LR™
Synonym : Disintegrant
HS Code : 3505 10 90
Shelf Life : 4 years

Technical Specification

Appearance : White to off-white free -flowing powder
Solubility : Readily soluble in water to form colloidal solution
pH (2% in water at 25°C) : 5.50 - 7.50
FTIR (KBr disc) : Matches with the standard pattern
Sodium chloride : <= 10.00%
Loss on drying (at 105°C) : <= 10.00%
Assay (sodium content, on dry basis) : 2.80 - 4.50%

Risk and Safety Information

Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods