

GRM7520 Sodium-DL-tartrate anhydrous**Product Number Packing**

GRM7520 : 500G

Product Information

Product Code : GRM7520
Product Name : Sodium-DL-tartrate anhydrous
Molecular Formula : $\text{Na}_2\text{C}_4\text{H}_4\text{O}_6$
Molecular Weight : 194.05
CAS No. : 868-18-8
HS Code : 2918 13 00
Shelf Life : 4 years

Technical Specification

Appearance : Colourless to white crystals or powder
Solubility : 100 mg soluble in 1 mL of water
Chloride (Cl) : $\leq 0.01\%$
Iron (Fe) : $\leq 0.004\%$
Lead (Pb) : $\leq 0.002\%$
Sulfate (SO_4) : $\leq 0.05\%$
Assay (NT) : 99.00 - 102.00%

Risk and Safety InformationStorage Temperature($^{\circ}\text{C}$) : Store below 30°C **Transport Information**

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods