
Strontium carbonate, Hi-LR™

GRM7523

Product Identifier

CAS No.	:	1633-05-2
EC No.	:	216-643-7
Molecular Formula	:	SrCO ₃
Molecular Weight	:	147.63
HS Code	:	2836 92 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white powder
Solubility	:	33.3 mg soluble in 1 mL of dilute acids
Chloride (Cl)	:	<= 0.005%
Iron (Fe)	:	<= 0.002%
Sulfate (SO ₄)	:	<= 0.001%
Assay (EDTA Titration)	:	min. 99.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	NWG