

Tetrabutylammonium hydrogensulfate, Hi-AR[®]

GRM7531


Product Identifier

CAS No.	:	32503-27-8
EC No.	:	251-068-5
Molecular Formula	:	C ₁₆ H ₃₇ NO ₄ S
Molecular Weight	:	339.53
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to off-white crystals or powder or lumps
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Assay (NaOH Titration)	:	98.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H314- H412
Precautionary Statement(s)	:	P260- P273- P280- P301+P330+P331- P305+P351+P338- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	2