

Sodium benzoate

GRM754


Product Identifier

CAS No.	:	532-32-1
EC No.	:	208-534-8
Molecular Formula	:	C ₇ H ₅ O ₂ Na
Molecular Weight	:	144.10
HS Code	:	2916 31 40
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or granules or powder or flakes
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Loss on drying (at 105°C, 1 hr)	:	<= 1.0%
Assay (NT, on dry basis)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P280- P305+P351+P338- P337 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH6650000
WGK	:	1