

---

## Gallotannin, Hi-LR®

**GRM7541**

### Product Identifier

---

CAS No.	:	1401-55-4
EC No.	:	215-753-2
Molecular Formula	:	C <sub>76</sub> H <sub>52</sub> O <sub>46</sub>
Molecular Weight	:	1701.20
Synonym	:	Tannic acid powder, Hi-LR®
HS Code	:	3201 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Off-white to light-brown powder or solid
Solubility	:	33.3 mg soluble in 1 mL of ethanol
Loss on drying (at 105°C, 2hr)	:	<= 12.0%
Residue on ignition	:	<= 1.0%
Sugars, dextrin	:	Passes test

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WW5075000
WGK	:	2