

---

## Tungsten (VI) oxide, Hi-LR™

**GRM7580**

---

### Product Identifier

CAS No.	:	1314-35-8
EC No.	:	215-231-4
Molecular Formula	:	WO <sub>3</sub>
Molecular Weight	:	231.84
HS Code	:	2815 20 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	Light yellow to dark green crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of hot 10% sodium hydroxide solution
Assay (Gravimetry)	:	min. 98.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	YO7760000
WGK	:	3