

GRM761 Strontium chloride hexahydrate, Hi-LR™**Product Number Packing**

GRM761 : 500G

Product Information

Product Code : GRM761
Product Name : Strontium chloride hexahydrate, Hi-LR™
Molecular Formula : $\text{SrCl}_2 \cdot 6\text{H}_2\text{O}$
Molecular Weight : 266.62
CAS No. : 10025-70-4
EC No. : 233-971-6
HS Code : 2827 39 85
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder or granules or chunks
Solubility : 33.3 mg soluble in 1 mL of water
pH (5% in water, at 25°C) : 5.00 - 7.00
Barium (Ba) : $\leq 0.1\%$
Calcium (Ca) : $\leq 0.1\%$
Iron (Fe) : $\leq 0.002\%$
Assay (EDTA Titration/AT) : 98.00 - 103.00%

GHS Safety Information

Hazard Statement(s) : H318
Precautionary Statement(s) : P280-P305+P351+P338
Signal Word : Danger

Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 41
S-Phrase(s) : 26-39
WGK : 1
RTECS : WK8575000
Storage Temperature(°C) : Store below 30°C

Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods



Product Information

IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein