

## Victoria blue B

**GRM7614**


### Product Identifier

CAS No.	:	2580-56-5
EC No.	:	219-943-6
C.I. No.	:	44045
Molecular Formula	:	C <sub>33</sub> H <sub>32</sub> ClN <sub>3</sub>
Molecular Weight	:	506.08
HS Code	:	3204 13 52
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Dark blue to dark brown glistening crystals or powder
Solubility	:	20 ppm in water yields clear solution
Absorption maxima	:	615 - 619 nm (20 ppm in water)
Loss on drying (at 60°C, 2 hr)	:	<= 5.0%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H319
Precautionary Statement(s)	:	P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3