
Xanthan gum

GRM7618

Product Identifier

CAS No.	:	11138-66-2
EC No.	:	234-394-2
HS Code	:	3913 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to tan free flowing powder
Solubility	:	33.3 mg soluble in 1 mL of water to give highly viscous solution
Loss on drying (at 105°C, 1 hr)	:	<= 15.00%
Pyruvic acid	:	>= 1.50%
Viscosity (1% Xanthan gum in 1% KCl)	:	>= 1200 cps (Spindle No. 3, 60 RPM at 25°C)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1