

GRM764 Soapstone**Product Number Packing**

GRM764 : 500G

Product Information

Product Code : GRM764
Product Name : Soapstone
Synonym : Talcum powder; Hydrous magnesium silicate; Talc, fine powder
Molecular Formula : $3\text{MgO} \cdot 4\text{SiO}_2 \cdot \text{H}_2\text{O}$
Molecular Weight : 379.27
CAS No. : 14807-96-6
EC No. : 238-877-9
HS Code : 2526 20 00
Shelf Life : 4 years

Technical Specification

Appearance : White to off-white powder
Solubility : Insoluble in water
Chloride (Cl) : $\leq 0.025\%$
Iron (Fe) : $\leq 0.001\%$
Acidity or alkalinity : Passes test
Acid soluble substance : $\leq 2.0\%$
Water soluble substance : $\leq 10 \text{ mg}$
Carbonates : No effervescence is produced
Organic compounds : The residue obtained in the test for loss on drying is not more than slightly yellow or grey
Loss on drying (at 180°C, 1hr) : $\leq 1.0\%$

GHS Safety Information

Hazard Statement(s) : 332
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 20-37
S-Phrase(s) : 36
WGK : nwg
RTECS : WW2710000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)





Product Information

Transport Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein