
Jenner's stain, Hi-Cert™

GRM767

Product Identifier

CAS No.	:	62851-42-7
HS Code	:	3204 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Green to brown or greenish black crystals or powder
Solubility	:	10 ppm in methanol yields dark-blue to very dark-blue clear solution
Absorption maxima 1	:	512 - 532 nm (10 ppm in methanol)
Absorption maxima 2	:	643 - 663 nm (10 ppm in methanol)
Specific Absorbance 1	:	>= 600 (10 ppm in methanol at 512 - 532 nm)
Specific Absorbance 2	:	>= 1100 (10 ppm in methanol at 643 - 663 nm)
Loss on drying (at 105°C, 1 hr)	:	<= 10.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1