

Atropine

GRM7851


Product Identifier

CAS No.	:	51-55-8
EC No.	:	200-104-8
Molecular Formula	:	C ₁₇ H ₂₃ NO ₃
Molecular Weight	:	289.37
HS Code	:	2939 79 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to off-white crystals or powder
Solubility	:	100 mg soluble in 1 mL of methylene chloride
FTIR	:	Matches with the standard pattern
Melting range	:	114 - 118°C
Assay (NT/HPLC)	:	99.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H300- H330
Precautionary Statement(s)	:	P260- P264- P284- P301 + P310- P304 + P340- P330
UN No.	:	1544
Class	:	6.1
Packing Group	:	II
RTECS	:	CK0700000
WGK	:	3