

6-Benzylaminopurine

GRM789


Product Identifier

CAS No.	:	1214-39-7
EC No.	:	214-927-5
Molecular Formula	:	C ₁₂ H ₁₁ N ₅
Molecular Weight	:	225.25
Synonym	:	6-BAP
HS Code	:	2933 59 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to light yellow crystals or powder
Solubility	:	50 mg soluble in 1 mL of glacial acetic acid gives clear to slightly hazy solution
FTIR	:	Matches with the standard pattern
Melting range	:	226 - 236°C
Assay (NT/HPLC)	:	98.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AU6252200
WGK	:	1

Other Information

Plant growth hormone