

2,2'-Bipyridine, Hi-AR®/ACS

GRM791


Product Identifier

CAS No.	:	366-18-7
EC No.	:	206-674-4
Molecular Formula	:	C ₁₀ H ₈ N ₂
Molecular Weight	:	156.18
Synonym	:	2,2'-Dipyridyl, Hi-AR®/ACS
HS Code	:	2933 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	100 mg soluble in 1 mL of ethanol
Melting range	:	69 - 73°C
Clarity and absorbance of solution	:	Passes test
Sensitivity to ferrous iron	:	Passes test
Assay (NT/GC)	:	99.50 - 101.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301- H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P301 + P310- P312- P403 + P233- P501
UN No.	:	2811
Class	:	6.1
Packing Group	:	III
RTECS	:	DW1750000
WGK	:	1