
Butylated hydroxyanisole

GRM796

Product Identifier

CAS No.	:	121-00-6
EC No.	:	204-442-7
Molecular Formula	:	C ₁₁ H ₁₆ O ₂
Molecular Weight	:	180.24
Synonym	:	3-tert-Butyl-4-hydroxyanisole; BHA
HS Code	:	2909 50 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow with tan cast or pink crystals or powder or waxy solid
Solubility	:	100 mg soluble in 1 mL of ethanol
FTIR	:	Matches with the standard pattern
Melting range	:	48 - 67°C
Assay (GC)	:	min. 98.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	SK1575000
WGK	:	2