

Butylated hydroxytoluene

GRM797


Product Identifier

CAS No.	:	128-37-0
EC No.	:	204-881-4
Molecular Formula	:	C ₁₅ H ₂₄ O
Molecular Weight	:	220.35
Synonym	:	BHT
HS Code	:	2907 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellowish-white crystals or powder
Solubility	:	100 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Assay (GC)	:	>= 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H410
Precautionary Statement(s)	:	P273- P391- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	GO7875000
WGK	:	2