

GRM8044**Cetyl palmitate (Mixture of C₁₄ - C₁₈ esters)****Product Number****Packing**

GRM8044 : 500G

Product Information

Product Code : GRM8044
Product Name : Cetyl palmitate (Mixture of C₁₄ - C₁₈ esters)
Synonym : Hexadecyl hexadecanoate; Palmitic acid palmityl ester
Molecular Formula : C₃₂H₆₄O₂
Molecular Weight : 480.85
CAS No. : 540-10-3
EC No. : 208-736-6
HS Code : 2915 70 00
Shelf Life : 4 years

Technical Specification

Appearance : White to off-white crystals or powder or leaflets
Solubility : 100 mg soluble in 1 mL of chloroform
Melting range : 46 - 56°C
Acid value : <= 4.00
Saponification value : 109 - 130
Iodine : <= 1.00
Assay (GC, mixture of esters) : min. 98.00%

Risk and Safety Information

WGK : nwg
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods