

GRM8051 4-Chloro-2-aminophenol**Product Number Packing**

GRM8051 : 500G

Product Information

Product Code : GRM8051
Product Name : 4-Chloro-2-aminophenol
Synonym : 5-Chloro-2-hydroxyaniline; 2-Amino-4-chlorophenol
Molecular Formula : C₆H₆NOCl
Molecular Weight : 143.57
CAS No. : 95-85-2
EC No. : 202-458-9
HS Code : 2922 29 00
Shelf Life : 4 years

Technical Specification


Appearance : White to brown crystals or powder or chunks
Solubility : 33.3 mg soluble in 1 mL of methanol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 134 - 142°C
Assay (GC/NT) : 99.00 - 103.00%

GHS Safety Information

Hazard Statement(s) : H302-H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22-36/37/38
S-Phrase(s) : 26-36
WGK : 3
RTECS : SJ5700000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 2673
Class : 6.1
Packaging Group : 2



Product Information

Marine Pollutant	:	No
ADR/RID	:	2673 6.1/PG 2
IMDG	:	2673 6.1/PG 2
IATA	:	2673 6.1/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein