

---

## Ehrlich's reagent

**GRM809**

### Product Identifier

---

CAS No.	:	100-10-7
EC No.	:	202-819-0
Molecular Formula	:	C <sub>9</sub> H <sub>11</sub> NO
Molecular Weight	:	149.19
Synonym	:	p-Dimethylaminobenzaldehyde
HS Code	:	2922 39 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to beige to green crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern.
Melting range	:	70 - 75°C
Assay (NT/GC)	:	99.00 - 103.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	CU5775000
WGK	:	1