
Florisil®

GRM817

Product Identifier

CAS No.	:	1343-88-0
EC No.	:	215-681-1
Molecular Formula	:	MgO ₃ Si
Molecular Weight	:	100.39
Synonym	:	Magnesium silicate activated, 60 - 100 mesh
HS Code	:	2839 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white fine powder
Solubility	:	Insoluble in water
Particle size 1	:	min. 98% passing through 60 mesh (250 µm)
Particle size 2	:	min. 70% passing through 100 mesh (149 µm)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	NWG