

## Indole

GRM824


### Product Identifier

CAS No.	:	120-72-9
EC No.	:	204-420-7
Molecular Formula	:	C <sub>8</sub> H <sub>7</sub> N
Molecular Weight	:	117.15
Synonym	:	2,3-Benzopyrrole
HS Code	:	2933 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to light pink crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Melting range	:	50 - 54°C
Assay (GC)	:	min. 99.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H311- H315- H318- H335- H400
Precautionary Statement(s)	:	P261- P273- P280- P305+P351+P338- P361 + P364- P403 + P233
UN No.	:	2811
Class	:	6.1
Packing Group	:	III
RTECS	:	NL2450000
WGK	:	1