

RM8245 N,N-Diisopropylethylamine**Product Number Packing**

RM8245	: 100ML
RM8245	: 500ML

Product Information

Product Code	: RM8245
Product Name	: N,N-Diisopropylethylamine
Synonym	: N-Ethyldiisopropylamine
Molecular Formula	: C ₈ H ₁₉ N
Molecular Weight	: 129.24
CAS No.	: 7087-68-5
EC No.	: 230-392-0
Shelf Life	: 4 years

Technical Specification

Appearance	: Colourless to yellowish liquid
Solubility	: 33.3 mg soluble in 1 mL of ethanol
FTIR (Liquid film)	: Matches with the standard pattern
Refractive index (n _{20/D})	: 1.412 - 1.415
Density (at 25°C)	: ~0.742 g/mL
Assay (GC)	: min. 98.00 %

GHS Safety Information

Hazard Statement(s)	: H225-H301-H314-H412
Precautionary Statement(s)	: P210-P273-P280-P301+P310-P305+P351+P338-P310
Signal Word	: Danger
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 11-22-34-52/53
S-Phrase(s)	: 16-26-36/37/39-45-60
WGK	: 2
Flash Point(°F)	: 50 °F
Flash Point(°C)	: 10 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

**Transport Information**

UN No.	: 2733
--------	--------



Product Information

Class	:	3
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	2733 3/PG 2
IMDG	:	2733 3/PG 2
IATA	:	2733 3/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein