
Kojic acid

GRM826

Product Identifier

CAS No.	:	501-30-4
EC No.	:	207-922-4
Molecular Formula	:	C ₆ H ₆ O ₄
Molecular Weight	:	142.11
HS Code	:	2932 99 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to beige-brown crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	148 - 157°C
Assay (NaOH Titration)	:	98.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UQ0875000
WGK	:	3