

2-Isopropoxyethanol

GRM8525


Product Identifier

CAS No.	:	109-59-1
EC No.	:	203-685-6
Molecular Formula	:	C ₅ H ₁₂ O ₂
Molecular Weight	:	104.15
HS Code	:	2909 44 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to light yellow liquid.
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Density (at 25°C)	:	0.903 - 0.904 g/mL
Assay (GC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H226- H312+H332- H315- H319
Precautionary Statement(s)	:	P210- P261- P280- P312- P370 + P378- P501
UN No.	:	3271
Class	:	3
Packing Group	:	III
RTECS	:	KL5075000
WGK	:	1

Other Information

Refractive index (n_{20/D}): 1.4050 to 1.4150