

Lead (II) iodide

GRM8547

Product Identifier

CAS No.	:	10101-63-0
EC No.	:	233-256-9
EC Index No	:	082-001-00-6
Molecular Formula	:	PbI ₂
Molecular Weight	:	461.01
HS Code	:	2827 60 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to orange powder.
Solubility	:	Insoluble in water
Assay (EDTA Titration)	:	98.50 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302+H332- H360Df- H373- H410
Precautionary Statement(s)	:	P201- P260- P280- P301+P312+P330- P308 + P313
UN No.	:	2291
Class	:	6.1
Packing Group	:	III
RTECS	:	Not available
WGK	:	3