

## 2-Methoxyethanol

**GRM8601**


### Product Identifier

CAS No.	:	109-86-4
EC No.	:	203-713-7
Molecular Formula	:	C <sub>3</sub> H <sub>8</sub> O <sub>2</sub>
Molecular Weight	:	76.09
HS Code	:	2909 44 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Colourless to faint yellow liquid
Solubility	:	1 mL miscible in 1 mL of water
FTIR	:	Matches with the standard pattern
Refractive index (n <sub>20/D</sub> )	:	1.400 - 1.404
Density (at 25°C)	:	0.960 - 0.970 g/mL
Assay (GC)	:	min. 99.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H226- H302+H312+H332- H360FD- H370- H373
Precautionary Statement(s)	:	P201- P210- P260- P280- P308 + P311- P370 + P378
UN No.	:	1188
Class	:	3
Packing Group	:	III
RTECS	:	KL5775000
WGK	:	3