

GRM8613 Methyl acetate

Product Number Packing

GRM8613 : 500ML

Product Information

Product Code : GRM8613
Product Name : Methyl acetate
Molecular Formula : $C_3H_6O_2$
Molecular Weight : 74.08
CAS No. : 79-20-9
Shelf Life : 4 years

Technical Specification

Appearance : Colourless clear liquid, having pleasant odour.
Solubility : Soluble in water to give clear to slightly hazy solution. Miscible with alcohol, ether.
Assay (GC) : min. 98.50 %

Risk and Safety Information

R-Phrase(s) : 11-36-66-67
S-Phrase(s) : 16-26-29-23
RTECS : AI9100000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

UN No. : 1231
Class : 3
Packaging Group : 2
ADR/RID : 1231 3/PG 2
IMDG : 1231 3/PG 2
IATA : 1231 3/PG 2