

Methyl acetoacetate

GRM8614


Product Identifier

CAS No.	:	105-45-3
EC No.	:	203-299-8
Molecular Formula	:	C ₅ H ₈ O ₃
Molecular Weight	:	116.12
HS Code	:	2918 30 40
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless liquid
Solubility	:	1 mL miscible in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n _{20/D})	:	1.414 - 1.424
Density (at 25°C)	:	1.071 - 1.081 g/mL
Assay (GC)	:	min. 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P280- P305+P351+P338- P337 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AK5775000
WGK	:	1