

Methylamine 37 - 43%, aqueous solution

GRM8616


Product Identifier

CAS No.	:	74-89-5
EC No.	:	200-820-0
Molecular Formula	:	CH ₅ N
Molecular Weight	:	31.06
HS Code	:	2921 11 90
Storage	:	Below 30°C
Shelf life	:	2 years

Technical Specification

Appearance	:	Colorless liquid
Solubility	:	1 mL miscible with 1 mL of water
Refractive index (n 20/D)	:	1.3650 - 1.3750
Density (at 25°C)	:	0.892 - 0.902 g/mL
Assay (HCl Titration)	:	37.0 - 43.0%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H224- H302+H332- H314- H335
Precautionary Statement(s)	:	P210- P260- P280- P303+P361+P353- P305+P351+P338- P370 + P378
UN No.	:	1235
Class	:	3(8)
Packing Group	:	II
RTECS	:	PF6300000
WGK	:	1