
Methyl oleate, Technical grade

GRM8651

Product Identifier

CAS No.	:	112-62-9
EC No.	:	203-992-5
Molecular Formula	:	C ₁₉ H ₃₆ O ₂
Molecular Weight	:	296.49
HS Code	:	2916 15 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow liquid
Solubility	:	1 mL miscible in 1 mL of anhydrous alcohol
Acid value	:	<= 5
Iodine value	:	85 - 120
Saponification value	:	180 - 195
Assay (GC)	:	min. 68.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RK0895000
WGK	:	1