
Glyoxal solution (40 wt.% in water)

GRM866


Product Identifier

CAS No.	:	107-22-2
Molecular Formula	:	C ₂ H ₂ O ₂
Molecular Weight	:	58.04
HS Code	:	2912 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow liquid (may darken in storage)
Solubility	:	1 mL miscible in 1 mL of water
Assay (NaOH Titration)	:	38.00 - 42.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H317- H319- H332- H341
Precautionary Statement(s)	:	P201- P280- P302 + P352- P304+P340+P312- P305+P351+P338- P308 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	MD2700000
WGK	:	3