

Sodium lauryl sulphate, Hi-LR[®]

GRM886


Product Identifier

CAS No.	:	151-21-3
EC No.	:	205-788-1
Molecular Formula	:	C ₁₂ H ₂₅ NaO ₄ S
Molecular Weight	:	288.38
Synonym	:	Sodium dodecyl sulphate, Hi-LR [®]
HS Code	:	2920 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water gives clear to opalescent solution
Assay (Hyamine Titration)	:	min. 85.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H228- H302- H311- H315- H318- H335
Precautionary Statement(s)	:	P210- P280- P305+P351+P338- P310- P370 + P378- P403 + P233
UN No.	:	1325
Class	:	4.1
Packing Group	:	III
RTECS	:	WT1050000
WGK	:	2