

---

## Sodium glycocholate hydrate

**GRM8907**

### Product Identifier

---

CAS No.	:	338950-81-5
EC No.	:	212-730-9
Molecular Formula	:	C <sub>26</sub> H <sub>42</sub> NO <sub>6</sub> Na.xH <sub>2</sub> O
Molecular Weight	:	487.60 (Anhydrous basis)
HS Code	:	2915 70 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to tan crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
Water (K.F.)	:	<= 10.00%
Assay (NT/HPLC, on anhydrous basis)	:	97.00 - 102.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3