

## m-Toluic acid

**GRM8997**

### Product Identifier

CAS No.	:	99-04-7
EC No.	:	202-723-9
Molecular Formula	:	C <sub>8</sub> H <sub>8</sub> O <sub>2</sub>
Molecular Weight	:	136.15
HS Code	:	2916 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to beige crystals or powder or flakes
Solubility	:	33.3 mg soluble in 1 mL of ethanol
FTIR	:	Matches with the standard pattern
Assay (NaOH Titration)	:	>= 98.00%

### Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H318
Precautionary Statement(s)	:	P280- P305+P351+P338- P310
UN No.	:	Not Dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XU1200000
WGK	:	1