



# Product Information

## GRM9354 Carboxymethyl cellulose calcium salt

---

| Product Number | Packing |
|----------------|---------|
|----------------|---------|

---

|         |        |
|---------|--------|
| GRM9354 | : 500G |
|---------|--------|

---

### Product Information

---

|              |  |
|--------------|--|
| Product Code | : GRM9354                              |
| Product Name | : Carboxymethyl cellulose calcium salt |
| CAS No.      | : 9050-04-8                            |
| Shelf Life   | : 4 years                              |

---

### Technical Specification

---

|                                 |                         |
|---------------------------------|-------------------------|
| Appearance                      | : White to cream powder |
| Solubility                      | : Insoluble in water    |
| Identification test             | : Passes test           |
| Loss on drying (at 105°C, 4 hr) | : <= 10.00%             |
| Residue on ignition             | : <= 20.00%             |

---

### Risk and Safety Information

---

|                         |                    |
|-------------------------|--------------------|
| Storage Temperature(°C) | : Store below 30°C |
|-------------------------|--------------------|

---

### Transport Information

---

|                  |                       |
|------------------|-----------------------|
| Marine Pollutant | : -                   |
| ADR/RID          | : Not Dangerous Goods |
| IMDG             | : Not Dangerous Goods |
| IATA             | : Not Dangerous Goods |

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein