
Leishman's stain, Practical grade (For Microscopy)

GRM949

Product Identifier

CAS No.	:	12627-53-1
EC No.	:	235-732-1
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Dark green to black powder
Solubility	:	20 ppm solution in methanol yields blue solution with green fluorescence
Absorption maxima 1	:	515 - 530nm (20 ppm in methanol)
Absorption maxima 2	:	640 - 650nm (20 ppm in methanol)
Loss on drying (at 105°C, 1 hr)	:	<= 8.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1