
May-Grunwald's stain, Practical grade

GRM953


Product Identifier

HS Code : 3212 90 90
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : Green to dark blackish green crystals or powder
Solubility : 20 ppm in methanol yields clear blue solution with greenish fluorescence
Absorption maxima 1 : 519 - 525 nm (20 ppm in methanol)
Absorption maxima 2 : 648 - 654 nm (20 ppm in methanol)
Loss on drying (at 105°C, 1 hr) : <= 10.00%

Safety Information

Hazard Pictogram(s) : 

Signal Word : Warning
Hazard Statement(s) : H302- H317- H319
Precautionary Statement(s) : P261- P280- P333 + P313- P337 + P313- P362 + P364- P501
UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : Not available
WGK : 1