
Potassium humate

GRM9781

Product Identifier

CAS No. : 68514-28-3
HS Code : 3808 93 40
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : Dark brown to black crystals or granules or flakes or powder
Solubility : 33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C) : 8.00 - 10.00%
Humic acid (Based on carbon content by T) : $\geq 65.00\%$

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : Not available
WGK : 1